

Thermal Analysis and Nano-Thermodynamics

Dr. Rakesh Kr Singh


Post-Doc, Ph.D

Thermal analysis

- **Thermal analysis** is a branch of **materials science** where the properties of materials are studied as they change with **temperature**. While engaging different techniques to inspect the measurement, both the types of instruments are employed to measure heat flow and this seems to be certified as DSC under the **ICTAC** (International Confederation for Thermal Analysis and Calorimetry) definition

Basic Principles of Thermal Analysis

- ▶ Modern instrumentation used for thermal analysis usually consists of the following parts:
 - ▶ sample holder/compartment for the sample
 - ▶ sensors to detect/measure a property of the sample and the temperature
 - ▶ an enclosure within which the experimental parameters (temperature, speed, environment) may be controlled
 - ▶ a computer to control data collection and processing



Thermal Analysis Techniques

- ▶ IUPAC definition - a group of techniques in which a physical property is measured as a function of temperature, while the sample is subjected to a controlled temperature programme (heating, cooling or isothermal).
- ▶ A range of techniques e.g.:
 - ▶ Differential Thermal Analysis (DTA) – temperature
 - ▶ Differential Scanning Calorimetry (DSC) – energy
 - ▶ Thermogravimetric Analysis (TGA) – mass
 - ▶ Thermomechanical Analysis (TMA) – dimensions

IUPAC: International Union of Pure and Applied Chemistry

Thermal Analysis

- Differential Scanning Calorimetry (DSC)
 - Measure heat absorbed or liberated during heating or cooling
- Differential Thermal Analysis (DTA)
 - They are use for thermal investigation where thermal change can be observed and characterised
- Thermal Gravimetric Analysis (TGA)
 - Measure change in weight during heating or cooling
- Thermomechanical Analysis (TMA)
 - Measure change in dimensions during heating or cooling


Thermogravimetric Analysis (TGA)- Basic Principle

TGA measures the amount and the rate of weight change of a material with respect to temperature or time in controlled environments. A TGA consists of three major parts a furnace-microgram balance, auto sampler and thermocouple.


- Nitrogen is the most common gas used to purge samples in TGA due to its inert nature.
- The furnace can raise the temperature as high as 1000°C which is made of quartz.
- The auto sampler helps to load the samples on to the microbalance.
- The thermocouple sits right above the sample.
- Care should be taken at all times that the thermocouple is not in touch with the sample which is in a platinum pan.

Thermogravimetric Analysis (TGA)


- ▶ A technique measuring the variation in mass of a sample undergoing temperature scanning in a controlled atmosphere
- ▶ Thermobalance allows for monitoring sample weight as a function of temperature
- ▶ The sample hangs from the balance inside the furnace and the balance is thermally isolated from the furnace



Thermogravimetric Analysis (TGA)



Examples of TGA Curves



TGA curves of crystalline and amorphous substance

Applications of TGA

- Thermal stability of materials
- Oxidative stability of materials
- Estimated lifetime of a product
- Decomposition Kinetics of materials
- The effect of reactive or corrosive atmosphere on materials
- Moisture and volatiles contents on materials.
- TGA reveals changes of a sample due to weight, whereas DTA and DSC reveal changes not related to the weight (mainly due to phase transitions)

Differential Thermal Analysis(DTA)

Sample holder: Sample and reference cells


Sensors: Thermocouples, one for the sample and one for the reference

Furnace: Block containing sample and reference cells

Temperature controller: Controls temperature program

Advantages:



- instruments can be used at very high temperatures
- instruments are highly sensitive
- flexibility in sample volume/form
- characteristic transition or reaction temperatures can be determined



Disadvantages:


- uncertainty of heats of fusion and transition temperatures

Differential Scanning Calorimetry (DSC)




- DSC measures differences in the amount of heat required to increase the temperature of a sample and a reference as a function of temperature

Sample containers and sampling




Differential Scanning Calorimeter




Principle Of DSC

Temperature,
Temperature
difference




Time


Typical DSC Curve




DSC Thermogram



Transitions in a DSC Curve




Interpretation of Thermal curves



- i. The sample undergoes no decomposition with loss of volatile products over the temperature range shown but solid phase transformation, melting, etc. can not be detected by TG,
- ii. The rapid initial mass loss is characteristic of desorption or drying. If it is true, then re-run the sample should result in type (i) curves,
- iii. Single stage decomposition,
- iv. Multi-stage decomposition with relatively stable intermediates: provide information on the temperature limit of stability of reactants and intermediate products and also stoichiometry,
- v. Multi-stage decomposition with no stable intermediate product. However heating-rate effect must be considered. At low heating rate, type (v) resemble type (iv). At high heating rate, type (iv) and (v) resemble type (iii) and lose all the details,
- vi. Gain in mass due to reaction with atmosphere, e.g. oxidation of metals,
- vii. Oxidation product decompose again at higher temperature; this is not often encountered.

DSC and Heating rate Formula



$$dH/dt = C_p * dT/dt + f(T,t)$$


Where $dH/dt \rightarrow$ heat flow measured by DSC

$C_p \rightarrow$ heat capacity or weight of the sample

$dT/dt \rightarrow$ heating rate

$f(T,t) \rightarrow$ time dependant or kinetic component

Typical Features of a DSC Trace






Thermodynamics of Materials at 450C

Dependence of magnetic and structural properties of Ni_{0.5} M_{0.5} Fe₂O₄ (M=Co, Cu) nanoparticles and Thermodynamics of prepared materials

Rakesh Kr Singh et al. *Int.J. of Engineering, Science and Technology*, Vol. 2, No. 8, 2010, pp. 73-79



- The properties of nanosize ferrite Magnetic materials crucially depend on the synthesis temperature. In citrate precursor method, 450C seems to be a critical temperature for growth of Ni-Co and Ni-Cu ferrites in pure phase. The growth of particles becomes much faster beyond this temperature. Crystalline structure changes from cubic spinel to tetragonal spinel as 450C temperature is crossed. The particle size, lattice constant and also magnetization parameters show sharp changes at 450C. This suggests that nucleation/growth/reaction mechanism(Thermodynamic) is different below 450C from that above it(Growth mechanism are shown in figure)



Thermodynamics of Materials at annealing temperature 450C

Rakesh Kr Singh, H.C. Verma et al, Ni-Zn Ferrite Nanoparticles and their Applications, manthan, Int. J.(2011)p

Single annealing temperature 450C for all our samples of nanocrystalline Ni-Zn mixed ferrite. The most interesting case seems to be with $\text{Ni}_{0.8}\text{Zn}_{0.2}\text{Fe}_2\text{O}_4$, where both the coercive field (116.10 Oe) and the saturation magnetization (52.18 emu/g) are largest. Values of saturation magnetization is higher than 50 emu/g have so far been achieved by using other methods, only through sintering at temperatures much above 450C. The Mössbauer spectrum shows that Fe occupies both the A and B sites in the sample and superparamagnetic fluctuations are not significant. So, this temperature is assumed as standard temperature.



Sample	H_c (Oe)	M_r (emu/g)	M_s (emu/g)	Squareness/Particle size (nm)
$\text{Ni}_{0.2}\text{Zn}_{0.8}\text{Fe}_2\text{O}_4$	22.42	0.5744	17.50	0.033/16
$\text{Ni}_{0.4}\text{Zn}_{0.6}\text{Fe}_2\text{O}_4$	1.53	0.1049	43.43	0.002/7
$\text{Ni}_{0.5}\text{Zn}_{0.5}\text{Fe}_2\text{O}_4$	90.86	5.658	40.54	0.140/16
$\text{Ni}_{0.6}\text{Zn}_{0.4}\text{Fe}_2\text{O}_4$	35.77	2.678	43.64	0.059/9
$\text{Ni}_{0.75}\text{Zn}_{0.25}\text{Fe}_2\text{O}_4$	93.10	4.610	38.94	0.118/9
$\text{Ni}_{0.8}\text{Zn}_{0.2}\text{Fe}_2\text{O}_4$	116.10	11.38	52.18	0.148/18

Sign up for a free, 7-day trial

[Publications home](#) [Journals](#) [Magazines](#) [Newspapers](#) [Reference works a](#)

Home » Publications » Reports, newsletters, and transcripts » Trade newsletters » Nanotechnology Weekly » February 2

[Save](#) [Export](#) [Print](#) [Cite](#)

New Findings from R.K. Singh and Co-Authors Describe Advances in Nanotechnology.



Newspaper

Nanotechnology Weekly

February 6, 2012 ([Hide copyright information](#))

COPYRIGHT 2009 NewsRX. This material is published under license from the publisher through the Gale Group, Farmington Hills, Michigan. All inquiries regarding rights or concerns about this content should be directed to [Customer Service](#).

 Like 0 Tweet [Permalink](#)

According to the authors of **recent research from Bihar, India**, “Two nano aluminate spinel materials (ZnAl_2O_4 and NiAl_2O_4) were synthesized by the Citrate precursor method. The citrate precursors consisting of co precipitated citrates of Zn^{2+} or Ni^{2+} and aluminum were first subjected to thermal analysis (TG-DSC) for determining the optimum temperature for annealing.”

“Two step decomposition was observed incorporating dehydration and formation of the aluminate. The second step gives an endo peak (-2937 J/g) at 356 A degrees C in the DSC curve of the co precipitated nickel (II) citrate-aluminum citrate gel in O_2 atmosphere. Kinetic/mechanistic analysis of the TG data has also been-----.

[To read the full text of this article and others like it, try us out for 7 days, FREE!](#)

High Beam Research

2002 establishments in the Chicago, United States
Discover specific articles and research from thousands of credible published sources.

Recently viewed items

Article: **New Findings from R.K. Singh, Bihar, India and Co-Authors Describe Advances in Nanotechnology...**

Thermal, XRD, and magnetization studies on ZnAl_2O_4 and NiAl_2O_4 spinels, synthesized by citrate precursor method and annealed at 450 and 650 °C

Rakesh K. Singh · A. Yadav · A. Narayan · Mukesh Chandra · R. K. Verma

 29th STAC-ICC Conference Special Chapter
© Akadémiai Kiadó, Budapest, Hungary 2011

Abstract Two aluminate spinel materials (ZnAl_2O_4 and NiAl_2O_4) were synthesized by the citrate precursor method. The citrate precursors consisting of coprecipitated citrates of Zn^{2+} or Ni^{2+} and aluminum were first subjected to thermal analysis (TG-DSC) for determining the optimum temperature for annealing. Two step decomposition was observed incorporating dehydration and formation of the aluminate. The second step gives an endo peak (-2937 J/g) at 356 °C in the DSC curve of the coprecipitated nickel(II) citrate-aluminum citrate gel in O_2 atmosphere. Kinetic/mechanistic analysis of the TG data has also been carried out and values of E_a , ΔS^\ddagger , ΔG^\ddagger , and A were approximated. On the basis of the findings, 450 °C has been chosen for annealing of the gels. Annealing has also been done at 650 °C for 1 h in muffle furnace in an attempt to obtain nanometric particles of aluminates (MAl_2O_4) {M = Ni, Zn} and to find out their magnetic properties which could render them useful for chemical sensing applications, etc. The TG-DSC curves of various powders which were obtained on annealing at the two temperatures did exhibit thermal instability when carried out in N_2 atmosphere. NiAl_2O_4 and ZnAl_2O_4 spinels (particle size 17 and 34 nm,

respectively) are obtained in pure crystalline phase at 650 °C. ZnAl_2O_4 prepared this way shows coercivity values of 470 and 58.37 G and NiAl_2O_4 , 107 and 23.24 G when annealed at 450 and 650 °C, respectively. ZnAl_2O_4 prepared by a polymer precursor method and annealed at 1000 °C, has earlier been reported to have coercivity value of 469 G. Thus, the citrate precursor method is good for the synthesis of ZnAl_2O_4 , producing single phase nanocrystalline powder of high quality and crystallinity. The value of magnetization was found to be small in the present case for the NiAl_2O_4 spinel obtained at 450 °C.

Keywords Nickel aluminate · Zinc aluminate · Annealing temperature · XRD pattern · Coercivity · Magnetization · Thermal stability · TG-DSC · Kinetic parameters · Nano particles · Spinel

Introduction


Aluminate content brings diverse applications to materials. Ceramic materials (cements, castable ceramics, bioceramics, and electroceramics) are normally formed from cubic crystal systems [1] which include garnets and spinels. Many rare earth aluminate garnets are used as laser host when doped with Nd(III) or Yb(III) and also as scintillators and phosphors. Upon doping, interesting properties have been observed. Eu^{2+} , R^{3+} -doping for example, brings thermoluminescence properties to calcium aluminate materials [2]. And $\text{BaMgAl}_{10}\text{O}_{17}:\text{Eu}^{2+}$ (BAM), is an important blue phosphor of plasma display panel (PDP) [3]. The rare earth aluminate glasses find use as alternatives to sapphire for use in infrared windows. Single crystals of lanthanum aluminate are known to have application as a substrate for deposition of thin films of the high temperature

At Present Position of **Dr. Rakesh Kumar Singh** (Post-Doc, Ph.D, M.Sc.)
Assistant Professor

Centre for Nanoscience & Nanotechnology (University Department)
Aryabhatta Knowledge University, Patna

Assistant Professor in Physics, **Patna Women's College**, Patna University (Aug. 2004-2013)

TG – DSC of the precursor prepared by Citrate method was carried out by NETZSCH STA 449 apparatus. The dehydration and curing was found to be completing in the range of 100-500°C. Thus the sintering temperature was chosen as 400°C for complete crystallization



Thermal Analysis Instrument Manufacturers-References

- Perkin Elmer Thermal Analysis Systems

<http://www.perkin-elmer.com/thermal/index.html>

- TA Instruments

<http://www.tainst.com/>

- Mettler Toledo Thermal Analysis Systems

<http://www.mt.com/>

- Rheometric Scientific

<http://www.rheosci.com/>

- Haake

<http://polysort.com/haake/>

- NETZSCH Instruments

<http://www.netzsch.com/ta/>

- SETARAM Instruments

<http://setaram.com/>

- Instrument Specialists, Inc.

<http://www.instrument-specialists.com/>

THANK YOU